



**NEW!**

# Nimbo

## COMPARE & SAVE!

A less expensive alternative to the Kaye® posture walker!



## Nimbo lightweight Posterior Safety Roller

- KA 1200 Tyke, Goldenrod Yellow, 1/cs
- KA 2200 Junior, Lime Green, 1/cs
- KA 3200 Youth, Cornflower Blue, 1/cs
- KA 4200 Young Adult, Midnight Blue, 1/cs

*The Wenzelite line of posterior safety rollers supports the user from behind requiring them to procure a more upright posture and it also promotes trunk extension.*

### FEATURES:

- Aluminum Frame in bright colors.
- Height adjustable.
- Soft rubber wheels adhere to any surface and will not allow the Nimbo to slide backwards.
- One directional rear wheels with aluminum ratchet and pin mechanism.
- Anti-Reverse override bracket disengages the one directional rear wheels and allows forward & reverse mobility.
- Front swivel wheels enhance maneuverability and ease of turning.
- Locking mechanism easily switches front wheels from swivel to non-swivel. Each wheel is independently adjusted.
- Skid resistant, non-marring and easy rolling wheels are ideal for indoor and outdoor use.
- Latex-free handgrips with flanged ends.
- Folds easily for convenient transporting.
- 3 year limited warranty.



Many physical therapists prefer posterior (reverse) walkers for children with Cerebral Palsy and other neurological, orthopedic and developmental disorders.

| ITEM #  | HANDLEBAR HEIGHT RANGE | OVERALL WIDTH | WIDTH INSIDE HANDLE GRIP | DEPTH OPENED AT BASE | PRODUCT WEIGHT | CARTON SHIPPING WEIGHT | WEIGHT CAPACITY |
|---------|------------------------|---------------|--------------------------|----------------------|----------------|------------------------|-----------------|
| KA 1200 | 14.5" - 18.5"          | 21"           | 13.5"                    | 25"                  | 9 LBS.         | 12 LBS.                | 100 LBS.        |
| KA 2200 | 19" - 25"              | 22"           | 13.5"                    | 29"                  | 11 LBS.        | 13 LBS.                | 100 LBS.        |
| KA 3200 | 23" - 30.5"            | 25"           | 15"                      | 33"                  | 12 LBS.        | 14 LBS.                | 200 LBS.        |
| KA 4200 | 28.5"-36"              | 25.5"         | 16"                      | 36"                  | 13 LBS.        | 15 LBS.                | 200 LBS.        |

mobility

posterior posture walker

nimbo

Wenzelite Re/hab

To order call toll free:

New York: 877-224-0946 • California: 800-371-2266

drive  
MEDICAL DESIGN & MANUFACTURING

## ACCESSORIES:

### Pelvic Stabilizer

KA 1050 for use with all Nimbo Posture Walkers, 1/bx

The pelvic stabilizer is adjustable in width and features back & lateral pads and a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and, therefore, cannot stay centered in the walker.



### Nimbo Forearm Platform

KA 1035 Forearm Platform, for use with all Nimbo Posture Walkers 1pr/bx

KA 1035 S Mounting Bracket, for use with all Nimbo Posture Walkers 1pr/bx

Constructed of lightweight aluminum, the Nimbo forearm platforms are height and depth adjustable. The handgrips extend forward or back to the position needed, supinate, pronate and are angle adjustable. Each side is individually adjustable and allows custom positioning for each user. The contoured armrest is molded in one piece with a flexible edge for added-- comfort and has Velcro® straps to prevent arms from slipping. They can be mounted to the mounting bracket either towards the inside or outside of the walker for lateral adjustment. Spacers are included for additional lateral range.



### Seat Harness

CE 1070 S for use with KA 1200, KA 2200 and KA 3200, 1/bx

CE 1070 L for use with KA 4200, 1/bx

Adjustable sling seat supports user that cannot stand completely on its own.

*\* For proper sizing, measure wrist to floor height.*



mobility  
safety roller accessories  
posterior posture

Wenzel  
lite Re/hab

To order call toll free:  
New York: 877-224-0946 • California: 800-371-2266

drive  
MEDICAL DESIGN & MANUFACTURING