

P.O. #: _____	Quote No: _____	Order Date: <u> </u> <u> </u> <u> </u> <u> </u> <u> </u> <u> </u>	Sales Rep: _____
Dealer: _____		Ship to: _____	
Address: _____		Address: _____	
City: _____	Prov: _____	Postal Code: _____	City: _____
			State.: _____
			Zip Code: _____
Tel: _____	Fax: _____	Tel: _____	Fax: _____
Tag: _____		Confirm my order via e-mail: <u> </u> <u> </u> <u> </u>	
Dealer Contact: _____		Send PDG e-News: <u> </u> <u> </u> <u> </u> E-mail: _____	

Bentley Frame Type - 20 Degrees of Tilt

50330 Bentley manual tilt wheelchair \$3995

Frame Options

- 50331 **Standard base** and seat frame NC
(Max. Weight Capacity 250 lbs)
- 50546 Transit Tie Down System (TTDS) ^{1,2} NC
- 50547 **Heavy duty package** (double gas strut) ^{3,4} \$300
(Max. weight capacity 450 lbs.)
- 50661 **Extended Base** c/w Heavy Duty Rear Stabilizers ⁵ \$600
- 51304 **Extended Base** c/w Rear Floor Jacks with triple contact pad ⁶ \$600

Client Information (Ensures optimal set-up of chair)

50548 Weight (lbs.): _____ Height (ft./in.): _____

Seat Width (Inches)

- 50549 14" seat width NC
- 50550 15" seat width NC
- 50551 16" seat width NC
- 50552 17" seat width NC
- 50553 18" seat width NC
- 50554 19" seat width NC
- 50555 20" seat width NC
- 50556 21" seat width \$195
- 50557 22" seat width \$195

Custom Seat Width (Requires 4 Week Delivery)

- 50558 23" seat width ³ \$600
- 50559 24" seat width ³ \$600
- 50560 25" seat width ³ \$600
- 50561 26" seat width ³ \$600
- 50562 27" seat width ³ \$600
- 50563 28" seat width ³ \$900
- 50564 29" seat width ³ \$900
- 50565 30" seat width ³ \$900
- 50566 31" seat width ³ \$900
- 50567 32" seat width ³ \$900

Seat Depth (Adjustable)

- 50568 16" seat depth NC
- 50569 17" seat depth NC
- 50570 18" seat depth NC
- 50571 19" seat depth NC
- 50572 20" seat depth NC
- 50573 21" seat depth NC
- 50574 22" seat depth NC

Frame Color

- 50575 Blue Shimmer STD
- 50576 Tracer Red NC
- 50577 Black Shimmer NC
- 50578 Burnt Orange NC
- 50579 Misty Purple NC
- 50580 Forest Green NC
- 50581 Custom color: \$300

Arm Type

- 50582 No armrests STD
- 50583 Dual post fixed height NC
- 50584 Dual post adjustable height \$295
8" - 11" 10" - 14" 12" - 16"
- 50587 3" armrest extension - forward ⁸ \$59
- 50588 3" armrest extension - rearward ⁸ \$59
- 50589 Locking flip-up armrests, ⁹ \$295
Set at 10" 11" 12" 13" 14"

Armrest Options

- 50594 Desk length armpad, 2" x 10" NC
- 50595 Full length armpad, 2" x 12" NC
- 50596 3" wide arm pads (pr) ¹⁰ - full desk \$49
- 50598 Gel armpad, 2" x 12" \$115
- 50599 Gel armpad, 3.5" x 12" \$115
- 50600 Contour arm trough - right (4.5" X 13") ¹⁰ \$85
- 50601 Contour arm trough - left (4.5" X 13") ¹⁰ \$85
- 50602 Flat hand pad - left (4.5" X 5") ¹¹ \$55
- 50603 Flat hand pad - right (4.5" X 5") ¹¹ \$55
- 50604 Clothing guards (black plastic) ⁸ \$205
- 50605 Armrest wideners (pr) ¹³ \$375

Seat Pan Options

- 50606 No seat pan STD
- 50607 Solid seat pan \$140

Seat Cushion Options

- 50608 No seat cushion STD
- 50609 MaTRx-Vi cushion: 14"w - 20"w ¹² \$590
Depth: 16" - 17" 18" - 19" 20" - 21"
- 50610 Spare MaTRx-Vi Cushion Cover: 14"w - 20"w \$115
- 50611 MaTRx-Vi Cushion: 21"w - 22"w ¹² \$645
Depth: 18" - 19" 20" - 21"
- 50612 Spare MaTRx-Vi Cushion Cover: 21"w - 22"w \$125

Footnotes

- 1** TTDS includes 4 factory installed transport brackets and a wheelchair -anchored lap belt. TTDS has been crash tested in accordance with ANSI/RESNA WC 19. See owner's manual for complete details.
- 2** TTDS not retrofittable to existing models. Not field serviceable. Not compatible with Attendant Wheel Locks or Heavy Duty Package.
- 3** Max weight capacity 450 lbs. Requires 20" composite or 24" heavy duty spoked wheel. Fixed axles only.
- 4** Order as "add on" for seat widths up to 22". Heavy duty package included on seat width 23" or wider.
- 5** Provides additional rearward stability for clients with excessive movement. Foot operated spring loaded release. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O₂ Holder.
- 6** Lifts rear wheels off ground for increased stability. Available with 12" and 20" rear wheels only, Triple contact pad standard. Not compatible with Attendant Wheel Locks, Rear Stabilizers, IV Pole, O₂ Holder.
- 7** Contact PDG for specification.
- 8** Not compatible with flip-up armrests.
- 9** Not compatible with clothing guards, back upholstery.
- 10** May interfere with rear wheels when tilted, depending on configuration.
- 11** Requires contour arm trough.
- 12** Cushion seat width supplied as per chair seat width ordered.

Dealer: _____ P.O.#: _____ Quote # _____

Caster Options (See chart below for seat height compatibility)

- 50615 6" x 1-1/4" urethane NC
- 50616 6" x 1-1/4" pneumatic \$110
- 51311 6" x 1-3/4" urethane \$75
- 50617 8" x 1-1/4" urethane¹⁴ NC
- 50618 8" x 1-1/4" pneumatic¹⁴ \$110
- 50619 8" x 2" urethane caster¹⁴ \$135
- 50620 8" x 2" pneumatic caster¹⁴ \$135

Rear Wheel Options (See chart below for seat height compatibility)

- 50621 12" composite wheel NC
- 50622 16" composite wheel NC
- 50623 20" composite wheel¹⁵ NC
- 51312 20" spoke wheel¹⁵ \$75
- 51461 22" composite wheel NC
- 51460 22" spoke wheel \$75
- 50624 24" composite wheel NC
- 51313 24" spoke wheel \$75
- 50625 24" Heavy Duty spoke wheel^{15, 16} \$350

Seat Height Options

(CIRCLE) desired seat-to-floor height in chart below.

All dimensions are approximate.

		Seat Height				
Rear Wheel	6" Casters					
	12"	13"	14"	15"	16"	NC
			17"	18"	\$60	
16"	15"	16"	17"	18"	NC	
			19"	20"	\$60	
20"	13"	14"	15"	16"	NC	
			17"	18"	\$60	
22"	14"	15"	16"	17"	NC	
			18"	19"	\$60	
24"	15"	16"	17"	18"	NC	
			19"	20"	\$60	
Rear Wheel	8" Casters					
16"	15"	16"	17"	18"	NC	
			19"	20"	\$60	
24"	15"	16"	17"	18"	NC	
			19"	20"	\$60	

Rear Tire Options

- 50637 Urethane tires STD
- 50638 Pneumatic rear tires NC
- 50639 Pneumatic rear tires with flat-free inserts \$200

Axle Type

- 50640 Fixed rear axles¹⁵ STD
- 50641 Quick release rear axles NC

Hand Rim Options

- 50642 Steel hand rims¹⁷ STD
- 50643 Plastic coated hand rims¹⁸ \$95
- 50644 Plastic coated projection hand rims¹⁹ \$225
- 50645 Select number of projections: 8__ 12__ 16__
- 50648 Select projection length: 3/4"__ 1"__
- 50650 Select projection angle: 0°__ 45°__ 90°__

Wheel Locks

- 50653 Push to lock wheel locks STD
- 50654 Pull to lock wheel locks NC
- 50655 Attendant only wheel lock (push to lock) \$85
- 50656 User wheel locks and attendant wheel locks¹⁷ \$150
- 50657 Wheel lock extensions (pair) \$60
- 50658 Hill holders \$350

Anti-Tip Options

- 50659 Rear anti-tippers²⁰ \$99
- 50660 Front anti-tippers \$300

Back Post Style

- 50662 Push Bar Back - 25" high STD
- 50663 Push Bar Back - 31" high \$35
- 50664 Cane Handle Back - 25" high^{21, 22} \$35
- 50665 Cane Handle Back - 31" high^{21, 22} \$35
- 50666 Adjustable height / angle push bar handle^{22, 23} \$395

Back Options

- 50667 Seat to back angle - set at: STD
90°__ 98°__ 106°__ 114°__ 122°__
- 50673 Backrest extension 10" high²⁴ \$150
- 50674 Fold down back²⁵ NC
- 50675 Spring back suspension²⁵ STD

Back Upholstery

- 50676 No back upholstery STD
- 50677 Standard back upholstery²⁶ \$125
Specify height: 21"__ 22"__ 23"__
- 50681 Short upholstery - 25" back frame only²⁶ \$125
Specify height: 18"__ 19"__ 20"__
- 50684 Tall upholstery - 31" back frame only²⁶ \$125
Specify height: 24"__ 25"__ 26"__
- 50688 Adjustable tension upholstery^{26, 31} \$250
Specify height: 21"__ 22"__ 23"__

MaTRx Back (www.motionconcepts.com)

- 50689 MaTRx PB: \$625
Width: 15"-17"__ 17"-19"__ 19"/20"__
MaTRx PB Height: 12"__ 16"__ 20"__
- 50695 MaTRx PB HD \$995
Width: 20"/21"__ 22"/23"__ 24"/25"__
Height: 16"__ 20"__
- 50698 MaTRx PB or PB HD spare cover \$115
- 50699 MaTRx PB deep: \$845
Width: 15"-17"__ 17"-19"__ 19"/20"__
Height: 12"__ 16"__ 20"__
- 50703 MaTRx PB deep HD \$1145
Width: 20"/21"__ 22"/23"__ 24"/25"__
Height: 16"__ 20"__
- 50702 MaTRx PB deep or PD deep HD spare cover \$125

Tilt / Recline Options

- 50706 Reclining backrest (manual) max. recline - 30°²⁷ \$1200

Headrest Options (www.stealthproducts.com)

- 50707 Headrest Mounting plate (bolt on)^{28, 29} \$75
Stealth Headrest Hardware Style
- 50708 Multi-axis hardware - HMO475³⁰ \$239
- 50709 Headrest Link Hardware -TWB480 hardware³⁰ \$306
Stealth Headrest Pad Size:
- 50710 6"__ 8"__ 10"__ 14"__ \$226

Footnotes

- ¹⁴ Compatible only with 16" and 24" rear wheels.
- ¹⁵ Required for Heavy Duty Package.
- ¹⁶ Not compatible with Pneumatic rear tires with Flat Free Inserts.
- ¹⁷ 20" or 24" rear wheels only.
- ¹⁸ Require a 4 to 6 week lead time.
- ¹⁹ Not compatible with 12" rear wheels.
- ²⁰ Required safety option on a manual tilt chair.
- ²¹ Seat widths 14" - 22" only.
- ²² Not compatible with headrest mounting plate.
- ²³ Not compatible with Cane Handle Back..
- ²⁴ Not compatible with adjustable height / angle push bar.
- ²⁵ Not compatible with reclining back.
- ²⁶ Back height measured from seat pan to top of upholstery.
- ²⁷ Not compatible with spring back suspension and fold down back.
- ²⁸ Not compatible with head rest extension or adjustable angle push bar.
- ²⁹ Required when mounting Stealth Headrest on Push Bar Back.
- ³⁰ Requires headrest mounting plate.

Dealer: _____ P.O.#: _____ Quote # _____

Attendant Tilt / User Operated Self - Tilt

- 50714 Attendant tilt only STD
- 50716 Self tilt only: right armrest __ left armrest __³¹ \$99
- 50718 Dual activator: backrest - right __ left __³¹ \$195
- 50720 Armrest - right __ left __³¹
- 50722 Trigger extension \$48

Front Rigging ³⁸

- 50723 Omit front rigging NC
- 50724 Sector blocks only (no front rigging supplied) NC
- 50725 Plug-in removable 70° front rigging STD
- 50726 Plug-in adjustable knee angle front rigging \$290
- 50727 Swing away removable 60° front rigging NC
- 50730 Swing away removable 70° front rigging NC
- 50731 Swing away elevating legrest \$355

Foot Plates

- 50472 Composite footplates STD
- 50733 Adjustable angle footplates \$140
- 50734 Oversize aluminum footplates³³ \$120
- 50735 One-piece **FIXED** footplate³⁴ \$350
- 50736 One-piece **FLIP-UP** footplate³⁵ \$395

Extension Tube

- 50737 Standard length (13" - 17") STD
- 50738 Long extension tubes (15" - 19") \$75
- 50739 Solid extension tubes - extra strength³⁵ \$99

Footrest/Legrest Accessories

- 50740 Adjustable length heel loops \$55
- 50741 **BodyPoint padded calf strap** ⁴¹ \$49
- 50742 Miller's dynamic footrest mounts³⁶ \$99

Positioning Straps (www.bodypoint.com)

- 50743 BodyPoint - center-pull auto style hip belt \$85
- 50744 BodyPoint - padded, center-pull auto style hip belt \$125
- 50745 BodyPoint - 4 point hip belt \$150
- 50746 BodyPoint - chest strap \$75

Miscellaneous Options

- 51314 **One arm drive** ³⁷ \$3,500
- 51315 **Protective lower padding** \$425
- 51316 **IV Pole** ³⁹ \$295
- 50747 O² holder ³⁹ \$295
- 50748 Vent tray c/w battery tray⁴⁰ call for quote

Custom Modifications

Contact PDG for any custom modification requirements.

Total Retail Value - This order \$ _____

Have all desired specs been checked? A complete spec sheet ensures your chair will arrive on time and as you expect. HCPCS codes supplied are for reference only and do not ensure coverage or payment.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product.

Footnotes

- 31** Requires headrest mounting plate.
- 33** Not compatible with solid or long extension tubes.
- 34** Available on plug-in front rigging only.
- 35** Not compatible with oversize footplates.
- 36** Not available with elevated legrests.
- 37** Restrictions apply. Extended delivery time. Contact PDG for more information.
- 38** On chair pricing applies for only 1 pair of front rigging. For additional front rigging prices, contact PDG Customer Service.
- 39** Not compatible with TTDS (Transit Tie Down System), Rear Stabilizers, Rear Floor Jacks, Attendant Wheel Locks.
- 40** Vent and battery model and dimension must be supplied. May restrict chair specifications.
- 41** Not available for sizes greater than 22" - sizes greater than 22" are not padded.