



Quickie IRIS

EIZ4, Effective November 2008

ACCOUNT #: _____ SHIP TO: _____

DATE: _____ NAME: _____

P.O. #: _____ ADDRESS: _____

BUYER: _____ CITY/ STATE: _____

MARKED FOR: _____ ZIP CODE: _____ PHONE NUMBER: _____

The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

BASE MODEL		
<input type="checkbox"/>	183WM56	Quickie IRIS with True Fit Growth \$1,999.00

CUSTOM FRAME COLOR		
<input type="checkbox"/>	183C108	Zebra (Choose base frame color) \$210.00 <i>Std with white unless other base color is chosen. NA with Kolorfusion base colors.</i>
<input type="checkbox"/>	183C41	Berry Vibrations \$210.00
<input type="checkbox"/>	183C20	Graphite \$210.00
<input type="checkbox"/>	183C42	Stars and Stripes \$210.00
<input type="checkbox"/>	183C43	Mossy Oak Camouflage \$210.00

PRE-SET TILT ANGLE (55° range) (Default in Bold)		
<input type="checkbox"/>	183TA5 E1161	-5 to 50° NC <i>Not available with heavy duty package and 12.5" STFH.</i>
<input type="checkbox"/>	183TA6 E1161	0 to 55° NC
<input type="checkbox"/>	183TA7 E1161	+5 to 60° NC <i>Not available with Heavy Duty Package.</i>

FRAME WIDTH		
<input type="checkbox"/>	183W5	14" \$95.00
<input type="checkbox"/>	183W6	15" \$95.00
<input type="checkbox"/>	183W7	16" NC
<input type="checkbox"/>	183W8	17" NC
<input type="checkbox"/>	183W9	18" Std
<input type="checkbox"/>	183W10	19" NC
<input type="checkbox"/>	183W11 E2201	20" \$300.00
<input type="checkbox"/>	183W12 E2201	21" \$300.00 <i>Not available with removable stroller handle extension.</i>
<input type="checkbox"/>	183W13 E2201	22" \$300.00 <i>Not available with removable stroller handle extension.</i>

PRE-SET TILT ANGLE (40° range) (Default in Bold)		
<input type="checkbox"/>	183TA4 K0009	-5 to 35° NC <i>NA with heavy duty package and 12.5" STFH.</i>
<input type="checkbox"/>	183TA1 K0009	0 to 40° NC
<input type="checkbox"/>	183TA10 E1161	+5 to 45° NC <i>Not available with Heavy Duty Package.</i>

FRAME OPTIONS		
<input type="checkbox"/>	183FS23	Standard Base and Seat Frame Std <i>Maximum Weight: 250 lbs.</i>
<input type="checkbox"/>	183TO1	Transit Option \$210.00 <i>Heavy Duty only available with 8" x 2" casters only</i>
<input type="checkbox"/>	183HV1	Heavy Duty Package (Max: 350 lbs) \$395.00
<input type="checkbox"/>	183FT46	Contracture Seat Frame NC <i>Not available with 14" or 15" frame depths. Not available with drop seat pan.</i>
<input type="checkbox"/>	183FS24 E1029	Vent/Battery Tray Base \$425.00 <i>NA with 14" and 15" frame width; minimum 17" STFH. Not available with Attendant Wheel Lock.</i>

FRAME DEPTH SET UP (CHOOSE OPTION A OR B)

A. FRAME DEPTH SET-UP BY QUICKIE

When selecting the Desired Seat Depth and Back Type, see **FIGURE 1.1** if additional charges apply to the Frame Depth. Default Measurements are in Bold.

Desired Seat Depth	Back Type					
	3	2	1	Flush	-1	-2
14"	NC	NC	NC	NC	N/A	N/A
15"	NC	NC	NC	NC	N/A	N/A
16"	NC	NC	NC	NC	NC	N/A
17"	\$300	NC	NC	NC	NC	NC
18"	\$300	\$300	NC	NC	NC	NC
19"	\$300	\$300	\$300	NC	NC	NC
20"	N/A	\$300	\$300	\$300	NC	NC
21"	N/A	N/A	\$300	\$300	\$300	NC
22"	N/A	N/A	N/A	\$300	\$300	\$300

FIGURE 1.1

TILT RELEASE OPTIONS		
<input type="checkbox"/>	183TT8	Foot Release - No Cables NC <i>NA with Vent and Battery Tray. NA with 21-22" depth and +5° set-ups.</i> <input type="checkbox"/> Left <input type="checkbox"/> Right (Facing from Rear)
<input type="checkbox"/>	183TT12	Backrest Dual trigger cable releases NC

FRAME COLOR		
<input type="checkbox"/>	183C2	Black NC <i>NA with Zebra.</i>
<input type="checkbox"/>	183C10	Black Opal NC
<input type="checkbox"/>	183C1	Blue NC
<input type="checkbox"/>	183C38	Blue Opal NC
<input type="checkbox"/>	183C12	Blue Velvet NC
<input type="checkbox"/>	183C36	Candy Blue NC
<input type="checkbox"/>	183C35	Candy Red NC
<input type="checkbox"/>	183C39	Candy Turquoise NC
<input type="checkbox"/>	183C6	Forest Green NC
<input type="checkbox"/>	183C4	Midnight Purple NC
<input type="checkbox"/>	183C37	Orange NC
<input type="checkbox"/>	183C5	Pearl Pink NC
<input type="checkbox"/>	183C40	Razzleberry NC
<input type="checkbox"/>	183C11	Toxic Green NC
<input type="checkbox"/>	183C8	Yellow NC

A. DESIRED SEAT DEPTH

See **FIGURE 1.1** and **FIGURE 1.2**, Dimension A. The seat surface length required to accommodate the client in a properly seated position.

<input type="checkbox"/>	183SD7	14" NC
<input type="checkbox"/>	183SD8	15" NC
<input type="checkbox"/>	183SD9	16" NC
<input type="checkbox"/>	183SD10	17" NC
<input type="checkbox"/>	183SD11	18" NC
<input type="checkbox"/>	183SD12	19" NC
<input type="checkbox"/>	183SD13	20" NC
<input type="checkbox"/>	183SD14	21" NC
<input type="checkbox"/>	183SD33	22" NC

A. USER'S SEATED HEIGHT (Default in Bold)

See **FIGURE 1.2**, Dimension **B**. Distance from the seat surface to the top of the head.

<input type="checkbox"/> 183SJ01	Less than 37"	NC
<input type="checkbox"/> 183SJ02	37" or Greater	NC




FIGURE 1.2


A. USER'S WEIGHT (Default in Bold)

<input type="checkbox"/> 183BO01	0 - 250 lbs	NC
<input type="checkbox"/> 183BO02	251 - 300 lbs	NC
<input type="checkbox"/> 183BO03	301 - 350 lbs	NC

A. BACK TYPE

See **FIGURE 1.1** and **FIGURE 1.3**, Dimension **C**. The location of the backrest surface in relation to the front of the backpost. The measurement is taken from the front surface of the backpost to the deepest area of the backrest in contact with the user's body.

<input type="checkbox"/> 183BQ01	+3"	NC
<input type="checkbox"/> 183BQ05	+2"	NC
<input type="checkbox"/> 183BQ02	+1"	NC
<input type="checkbox"/> 183BQ03	Flush	NC
<input type="checkbox"/> 183BQ04	-1"	NC
<input type="checkbox"/> 183BQ06	-2"	NC



Back Type	
1	3" forward (+3')
2	2" forward (+2')
3	1" forward (+1')
4	Flush
5	1" rearward (-1")
6	2" rearward (-2")

FIGURE 1.3

A. TOP OF CUSHION TO SEAT RAIL (Default in Bold)

See **FIGURE 1.4** Dimensions **A** & **B**. This measurement is taken from the top of the cushion to the upper surface of the seat rail. If the cushion has a wedged shape and/or has significant lateral pelvic contour, then determine the height of the cushion at the point where the ischial tuberosities would contact (not immerse) the cushion surface.

<input type="checkbox"/> 183CJ01	2 1/2" or Greater	NC
<input type="checkbox"/> 183CJ02	Less than 2 1/2"	NC
<input type="checkbox"/> 183CJ03	Below Seat Rail	NC




FIGURE 1.4

B. FRAME DEPTH RANGE

<input type="checkbox"/> 183FL80	Short Seat Frame 14" - 18"	NC
<input type="checkbox"/> 183FL81	Medium Seat Frame 16" - 20"	NC
<input type="checkbox"/> 183FL82	Long Seat Frame 18" - 22"	NC

B. FRAME DEPTH

<input type="checkbox"/> 183FD5	14"	NC
	<i>Recommended Short Seat Frame.</i>	
	<i>NA with Contracture Footrest System.</i>	
<input type="checkbox"/> 183FD6	15"	NC
	<i>Recommend Short Seat Frame.</i>	
	<i>NA with Contracture Footrest System.</i>	
<input type="checkbox"/> 183FD7	16"	Std
	<i>Recommend Short Seat Frame.</i>	
<input type="checkbox"/> 183FD8	17"	NC
	<i>Recommend Short Seat Frame.</i>	
<input type="checkbox"/> 183FD9	18"	NC
	<i>Recommend Medium Seat Frame.</i>	
<input type="checkbox"/> 183FD10	19"	NC
	<i>Recommend Medium Seat Frame.</i>	
<input type="checkbox"/> 183FD11	E2203 20"	\$300.00
	<i>Recommend Long Seat Frame</i>	
<input type="checkbox"/> 183FD16	E2203 21"	\$300.00
	<i>Recommend Long Seat Frame.</i>	
<input type="checkbox"/> 183FD17	E2204 22"	\$300.00
	<i>Recommend Long Seat Frame.</i>	

SEAT TYPE

<input type="checkbox"/> 183ST35	K0108 Standard Seat Pan (Growable)	\$95.00
	<i>For Growth Ranges please see FIGURE 1.4 below.</i>	
<input type="checkbox"/> 183ST36	K0108 Drop Seat Pan	\$95.00
	<i>NA with transit option or HD Package</i>	
	<i>Not available with Contracture Footrest System.</i>	
	<i>Only available with 14, 16, 18, 20, 22 inch seat widths</i>	

Frame Width	Seat Pan Growth
14"	14" - 16"
15"	15" - 16"
16"	16" - 18"
17"	17" - 18"
18"	18" - 20"
19"	19" - 20"
20"	20" - 22"
21"	21" - 22"
22"	NA

FIGURE 1.4

SEAT CUSHION

To order any of the following JAY cushions on-chair, please attach the on-chair cushion t

<input type="checkbox"/> 183JC01	JAY Basic	See OF
<input type="checkbox"/> 183JC02	JAY Easy	See OF
<input type="checkbox"/> 183JC03	JAY Soft Combi P	See OF
<input type="checkbox"/> 183JC04	JAY J2	See OF
<input type="checkbox"/> 183JC05	JAY J2 Deep Contour	See OF
<input type="checkbox"/> 183JC06	JAY Lite	See OF

BACKPOST (SELECT OPTION 1, 2, or 3)

1. STROLLER HANDLE BACK POSTS

<input type="checkbox"/>	183BP9	Low (18")	NC
<input type="checkbox"/>	183BP10	Medium (21")	Std
<input type="checkbox"/>	183BP11	Tall (24")	NC
<input type="checkbox"/>	183BP13	Removable Extension	\$85.00
<i>NA with reclining backrest or 21" and 22" frame widths.</i>			
<i>NA with Height adj. back with adjustable stroller handle.</i>			

2. HEIGHT ADJ. BACK WITH ADJ. STROLLER HANDLES

<input type="checkbox"/>	183PH8	Standard Push Handle	NC
<input type="checkbox"/>	183PH9	Long Push Handle	NC
<i>Recommended in combo w/ Vent and Battery Tray Option.</i>			
<input type="checkbox"/>	183B200	Standard Height Back cane 15" - 21"	NC
<input type="checkbox"/>	183B201	Tall Height Back cane 18" - 24"	NC

3. RECLINING BACKREST

Only available with 40 degree tilt range

<input type="checkbox"/>	183B5	E1225 Medium	\$850.00
<input type="checkbox"/>	183B9	E1225 Tall	\$850.00

JAY J3 BACK

Please see the J3 Back order form for part #'s.

<input type="checkbox"/>	183JB01	Jay J3 Back	See OF
<i>NA with O2 Holder.</i>			

ARMREST

<input type="checkbox"/>	183AA63	E0973 Dual Post Height Adjustable	\$230.00
<input type="checkbox"/>	183AA43	Dual Post Flip Back	\$95.00
<input type="checkbox"/>	183AA62	Dual Post	NC
<input type="checkbox"/>	183AA48	E0973 Dual Post Flip Back Height Adjust.	\$245.00
<input type="checkbox"/>	183AA7	E0973 Length Adjustable Locking Flip-Up	\$235.00
<i>Not available with heavy duty package.</i>			
<input type="checkbox"/>	183AA8	E0973 Height Adjustable Adult	\$230.00
<i>Not available with heavy duty package.</i>			
<input type="checkbox"/>	183AA9	E0973 Height Adjustable Kids	\$230.00
<i>Not available with heavy duty package.</i>			
<input type="checkbox"/>	183AA1	E0973 Adjustable Locking Flip-Up	\$230.00
<i>NA with heavy duty package or arm pads.</i>			
<input type="checkbox"/>	183AA14	Omit Armrest	NC

ARMPAD

<input type="checkbox"/>	183AM1	Desk Length	NC
<input type="checkbox"/>	183AM2	Full Length	NC
<input type="checkbox"/>	183AM9	Foam Short	NC
<i>Option for Length Adj Locking Flip Up Armrest.</i>			
<input type="checkbox"/>	183AM10	Foam Long	NC
<i>Option for Length Adj Locking Flip Up Armrest.</i>			

FOOTREST (CHOOSE OPTION 1,2 or 3)

1. HANGER MOUNTED FOOTREST

<input type="checkbox"/>	183H7	70° Swing-Away Front Mount	NC
<input type="checkbox"/>	183H8	80° Swing-Away Front Mount	NC
<input type="checkbox"/>	183H6	90° Swing-Away Front Mount	NC
<input type="checkbox"/>	183H9	70° Heavy-Duty Front Mount Lift-Off	\$95.00
<input type="checkbox"/>	183H10	80° Heavy-Duty Front Mount Lift-Off	\$95.00
<i>Not available with 4" casters.</i>			
<input type="checkbox"/>	183H11	90° Heavy-Duty Front Mount Lift-Off	\$95.00
<i>Not available with 4" casters.</i>			
<input type="checkbox"/>	183H23	E0990 90° Elevating	\$185.00
<i>Not available with 6" casters or 90° footboard.</i>			

HANGER MOUNTED FOOTREST EXTENSION TUBES

NA with 12.5" seat height

<input type="checkbox"/>	183E23	2"	\$30.00
<i>NA with 90° swing-away and lift-off hangers.</i>			
<input type="checkbox"/>	183E24	4"	\$30.00
<i>NA w/ 90° SA and lift-off hangers.</i>			

HANGER MOUNTED FOOTREST FOOTPLATES

<input type="checkbox"/>	183F4	K0040 K0038 Kids Angle Adjust. w/Leg Strap	\$150.00
<i>Available with 14" to 16" frame width only.</i>			
<input type="checkbox"/>	183F3	K0040 E0951 Adult Angle Adjust. w/Heel Loops	\$150.00
<input type="checkbox"/>	183F62	90° Footboard w/Leg Strap	\$150.00
<i>Available with 14" to 16" frame width only.</i>			
<input type="checkbox"/>	183A7	Omit Footplate	NC

2. EXTENSION TUBE MOUNTED FOOTRESTS

Footrests Mount to Extension Tubes

<input type="checkbox"/>	183H1	60° Swing-Away	NC
<input type="checkbox"/>	183H2	70° Swing-Away	STD
<input type="checkbox"/>	183H3	70°V Swing-Away	NC
<i>NA with 12.5" seat to floor height.</i>			
<i>Available with 16" to 20" frame width only.</i>			
<i>Offered with composite or angle adjustable footplate.</i>			
<i>With composite footplates NA with 8" casters.</i>			
<input type="checkbox"/>	183H14	E0990 Elevating	\$185.00
<i>NA w/ locking angle adj. footplate</i>			
<input type="checkbox"/>	183H24	K0053 Articulating Adult	\$275.00
<i>NA w/ locking angle adj. footplate or short ext. tubes.</i>			

EXTENSION TUBE MOUNTED FOOTRESTS EXTENSION TUBES

<input type="checkbox"/>	183E1	Short (11.5" - 14")	NC
<input type="checkbox"/>	183E2	Medium (14" - 17")	NC
<i>NA with 12.5" seat height</i>			
<input type="checkbox"/>	183E5	Long (17" - 19")	NC
<i>NA with 12.5" seat height</i>			
<input type="checkbox"/>	183E16	High Mount	\$55.00
<i>Only available with adult or kids angle adj. footplates.</i>			

EXTENSION TUBE MOUNTED FOOTREST FOOTPLATES

<input type="checkbox"/>	183F1	Composite w/Heel Loops	STD
<i>With 70° SA Ext. mount hangers NA with 8" casters.</i>			
<input type="checkbox"/>	183F5	Aluminum	NC
<i>Only available with elevating legrests.</i>			
<input type="checkbox"/>	183F3	K0040 E0951 Adult Angle Adjust. w/Heel Loops	\$150.00
<input type="checkbox"/>	183F61	K0040 Locking Angle Adjustable	\$165.00
<i>Not available with V hangers.</i>			
<i>Only available with 15-22" seat widths.</i>			
<input type="checkbox"/>	183F48	Platform Flip-Up w/Leg Strap	\$155.00
<i>Not available with short extension tubes.</i>			
<i>Only available with 60°, 70°, or 70° V hangers.</i>			
<input type="checkbox"/>	183F7	Omit Footplate and Ext. Tube	NC

FOOTREST OPTIONS

<input type="checkbox"/>	183IG2	Impact Guard Neoprene	\$85.00
<i>Only available w/ 60°, 70°, or 70° V ext. mount hangers.</i>			
<input type="checkbox"/>	183A110	E0952 Toe Loop (pair)	\$20.00
<i>Only available with adult angle adjustable footplates.</i>			
<input type="checkbox"/>	183H32	Omit Hanger	-\$15.00
<input type="checkbox"/>	183HL4	Omit Heel Loop	NC
<input type="checkbox"/>	183HL6	Omit Single Leg Strap	NC

FRONT SEAT HEIGHT MATRIX

PLEASE SEE OPTION LISTING BELOW FOR IMPORTANT NO-GO INFORMATION.

		CASTER SIZE				
		<input type="checkbox"/> 4" Semi-Pneumatic	<input type="checkbox"/> 5" x 1" Polyurethane <input type="checkbox"/> 5" x 1.5" Semi-Pneumatic	<input type="checkbox"/> 6" x 1" Polyurethane <input type="checkbox"/> 6" x 1.5" Semi-Pneumatic	<input type="checkbox"/> 8" x 1" Polyurethane <input type="checkbox"/> 8" x 1.5" Semi-Pneumatic	<input type="checkbox"/> 8" x 2" Pneumatic <input type="checkbox"/> 8" x 2" Pneumatic w/ Airless Insert
REAR WHEEL SIZE	<input type="checkbox"/> 12" Mag Wheel	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13.5" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0" <input type="checkbox"/> 17.5"	<input type="checkbox"/> 16.0" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0"	<input type="checkbox"/> 17" <input type="checkbox"/> 17.5" <input type="checkbox"/> 18"	<input type="checkbox"/> 17" <input type="checkbox"/> 18"
	<input type="checkbox"/> 12" Low Poly Wheel		<input type="checkbox"/> 16.0" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0" <input type="checkbox"/> 17.5"	<input type="checkbox"/> 16.0" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0"	<input type="checkbox"/> 16" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17"	<input type="checkbox"/> 16.5" <input type="checkbox"/> 17.5"
	<input type="checkbox"/> 16" Mag Wheel		<input type="checkbox"/> 13.5" <input type="checkbox"/> 14.0" <input type="checkbox"/> 14.5" <input type="checkbox"/> 15.0" <input type="checkbox"/> 18.0" <input type="checkbox"/> 18.5"	<input type="checkbox"/> 14.5" <input type="checkbox"/> 15" <input type="checkbox"/> 18" <input type="checkbox"/> 18.5" <input type="checkbox"/> 19"	<input type="checkbox"/> 18" <input type="checkbox"/> 18.5" <input type="checkbox"/> 19"	<input type="checkbox"/> 18" <input type="checkbox"/> 19"
	<input type="checkbox"/> 20" Mag Wheel		<input type="checkbox"/> 16.5" <input type="checkbox"/> 17.0" <input type="checkbox"/> 17.5"	<input type="checkbox"/> 16" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17"	<input type="checkbox"/> 16" <input type="checkbox"/> 16.5" <input type="checkbox"/> 17"	<input type="checkbox"/> 16.5" <input type="checkbox"/> 17"
	<input type="checkbox"/> 22" Mag Wheel		<input type="checkbox"/> 17.0" <input type="checkbox"/> 17.5" <input type="checkbox"/> 18.0" <input type="checkbox"/> 18.5"	<input type="checkbox"/> 17" <input type="checkbox"/> 17.5" <input type="checkbox"/> 18"	<input type="checkbox"/> 17" <input type="checkbox"/> 17.5" <input type="checkbox"/> 18"	<input type="checkbox"/> 17" <input type="checkbox"/> 18"
	<input type="checkbox"/> 24" Mag Wheel <input type="checkbox"/> 24" Spoke Wheel		<input type="checkbox"/> 18.0" <input type="checkbox"/> 18.5"	<input type="checkbox"/> 18" <input type="checkbox"/> 18.5" <input type="checkbox"/> 19"	<input type="checkbox"/> 18" <input type="checkbox"/> 18.5" <input type="checkbox"/> 19" <input type="checkbox"/> 19.5"	<input type="checkbox"/> 18" <input type="checkbox"/> 19"

3. CONTRACTURE FOOTREST SYSTEM

Multi-Adjustable Padded Hangers & Angle Adjustable Aluminum Footrests

- 183H26 Contracture Footrest System \$325.00
*Not available with drop seat pan.
Not available with 15" frame depth.*
- 183FP02 Outside Mount NC
- 183FP01 Inside Mount NC
- 183E25 Short

Select Extension Height Below

- 183EP1 10.5" Extension Height NC
- 183EP2 11.0" Extension Height NC
- 183EP3 11.5" Extension Height NC
- 183EP4 12.0" Extension Height NC
- 183EP5 12.5" Extension Height NC
- 183EP6 13.0" Extension Height NC
- 183E26 Regular

Select Extension Height Below

- 183EP5 12.5" Extension Height NC
- 183EP6 13.0" Extension Height NC
- 183EP7 13.5" Extension Height NC
- 183EP8 14.0" Extension Height NC
- 183EP9 14.5" Extension Height NC
- 183EP10 15.0" Extension Height NC

CASTER WHEELS

For Reference Only, Please check box on matrix to order

- 183CW44 K0072 4" Semi-Pneumatic NC
Not available with Heavy Duty Option
- 183CW53 E2222 5" x 1" Polyurethane NC
- 183CW7 K0072 5" x 1.5" Semi-Pneumatic \$85.00
- 183CW45 E2222 6" x 1" Polyurethane Std
- 183CW35 K0072 6" x 1.5" Semi-Pneumatic \$85.00
- 183CW46 E2222 8" x 1" Polyurethane NC
- 183CW39 K0072 8" x 1.5" Semi-Pneumatic \$85.00
- 183CW36 K0071 8" x 2" Pneumatic \$145.00
- 183CW37 E2214 8" x 2" Pneu w/Airless Insert \$225.00

REAR WHEEL

For Reference Only, Please check box on matrix to order

- 183RW13 12" Mag Wheel NC
*NA with hand rims.
With attendant wheel lock n/a with +5 - 60° tilt angle.
Only avail. w/ pneu, pneu w/ airless inserts or low poly tires.*
- 183RW14 16" Mag Wheel NC
*NA with hand rims.
OA with pneumatic or pneumatic w/ airless inserts.*
- 183RW3 20" Mag Wheel NC
NA with Kevlar tires or projection handrims.
- 183RW8 22" Mag Wheel NC
NA with Kevlar tires or projection handrims.
- 183RW1 24" Mag Wheel NC
- 183RW2 24" Spoke Wheel \$200.00

REAR WHEEL TIRE

- 183RT3 Full Polyurethane NC
- 183RT4 Low Polyurethane NC
Only available with 12" Mag Wheel
- 183RT1 E2211 Pneumatic NC
- 183RT2 E2213 Pneumatic Airless Insert \$70.00
- 183RT8 E2211 Kevlar \$95.00
Only available with 24" rims.

HANDRIM

12" and 16" wheels not available with hand rims.

- 183HR1 Aluminum Anodized NC
- 183HR4 E2205 Plastic Coated \$85.00
- 183HR7 Omit Handrim -\$20.00

PROJECTIONS

Not available on plastic coated. Only available with 24" rims.

- 183HR8 E0967 8 Oblique \$195.00
- 183HR10 E0967 12 Oblique \$195.00
- 183HR9 E0967 8 Vertical \$195.00
- 183HR11 E0967 12 Vertical \$195.00

AXLE PLATE

<input type="checkbox"/>	183AL33	Standard Axle Plate	NC
<input type="checkbox"/>	183AL17 E0958	One Arm Drive Axle Plate <i>NA with 12" mag, 16" mag and 24" spoke wheels.</i>	\$795.00
<input type="checkbox"/>	183OA2	One Arm Drive Standard Handrim	NC
<input type="checkbox"/>	183HM1	Right Hand Mount	NC
<input type="checkbox"/>	183HM2	Left Hand Mount	NC
<input type="checkbox"/>	183OA1	One Arm Drive Small Handrim <i>Not available with 20" mag wheels.</i>	NC
<input type="checkbox"/>	183HM1	Right Hand Mount	NC
<input type="checkbox"/>	183HM2	Left Hand Mount	NC

AXLE

<input type="checkbox"/>	183AX1	Quick Release Axles <i>Not available with one arm drive.</i>	NC
<input type="checkbox"/>	183AX17	Threaded Axles <i>NA with 24 "spoke wheels</i>	NC

WHEEL LOCK

<input type="checkbox"/>	183WL8	Push to Lock	NC
<input type="checkbox"/>	183WL25	Attendant Wheel Lock (Foot Lock) <i>NA with fold down extension handles. Not available with 1 arm drive</i>	\$50.00
<input type="checkbox"/>	183WL10	Omit Wheel Lock	NC

EXTENSION HANDLE

<input type="checkbox"/>	183EH5 E0961	Fold Down <i>Not available with attendant wheel lock.</i>	\$60.00
--------------------------	--------------	--	---------

ANTI-TIP

Included with vent tray option. Not ordered as separate line item.

<input type="checkbox"/>	183AT1 E0971	Rear Anti-Tip	\$135.00
--------------------------	--------------	---------------	----------

SEAT POUCH

<input type="checkbox"/>	183A8	Black	\$45.00
--------------------------	-------	-------	---------

BACKPACK

<input type="checkbox"/>	183A1	Black Utility Backpack	\$55.00
--------------------------	-------	------------------------	---------

EXTRA LEG STRAP

<input type="checkbox"/>	183A54	Black	\$30.00
--------------------------	--------	-------	---------

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	183A82	Caddy <i>Only available with extention tube mounted footrests.</i>	\$45.00
<input type="checkbox"/>	183A71 E2208	O2 Holder <i>NA w/ Height Adj Stroller Handles. NA w/ Reclining Backrest, or Transit Option.</i>	\$150.00
<input type="checkbox"/>	183A202 K0105	IV Pole <i>NA w/ transit option.</i>	\$125.00
<input type="checkbox"/>	183A208	6" footrest extensions <i>May interfere with specific configurations</i>	\$30.00

POSITIONING BELT

<input type="checkbox"/>	183A158 E0978	2" Dual Pull Squeeze Pad w/Clamp	\$85.00
<input type="checkbox"/>	183A141 E0978	1.5" Auto Buckle Universal w/Clamp	\$60.00
<input type="checkbox"/>	183A117 E0978	2" Aircraft Buckle w/Clamps	\$85.00
<input type="checkbox"/>	183A118 E0978	2" Aircraft Buckle Pad w/Clamp	\$115.00
<input type="checkbox"/>	183A120 E0978	Positioning Belt Clamps Only	\$35.00

SPOKE GUARDS

Not available with 12" and 16" mag wheels.

<input type="checkbox"/>	183A24 K0065	Clear	\$90.00
<input type="checkbox"/>	183A23 K0065	Black	\$90.00

Not available with projection hand rims.

TOUCH-UP PAINT

<input type="checkbox"/>	183A27	Frame	\$25.00
--------------------------	--------	-------	---------

WHEELCHAIR TRAY TABLE

Only available if armrests are included.

<input type="checkbox"/>	183A30 E0950	Small (13"-14")	\$200.00
<input type="checkbox"/>	183A31 E0950	Medium (15"-17")	\$200.00
<input type="checkbox"/>	183A32 E0950	Large (18"-20")	\$200.00
<input type="checkbox"/>	183A33 E0950	X-Large (21"-22")	\$200.00

SIDE GUARD

Only available with length adjustable locking flip-up and adjustable locking flip-up armrest styles.

<input type="checkbox"/>	183SG3	Plastic Regular	\$120.00
<input type="checkbox"/>	183SG4	Plastic Kids	\$120.00