



# Luminator Posterior Gait Trainer

- LT 3000 Adult Luminator Posterior Gait Trainer, 1/cs
- LT 3100 Pediatric Luminator Posterior Gait Trainer, 1/cs
- LT 3200 Tyke Luminator Posterior Gait Trainer, 1/cs

# Luminator Anterior Gait Trainer

- LT 2000 Adult Luminator Anterior Gait Trainer, 1/cs
- LT 2100 Pediatric Luminator Anterior Gait Trainer, 1/cs
- LT 2200 Tyke Luminator Anterior Gait Trainer, 1/cs



Posterior



Anterior

## FEATURES:

- Seat Harness – supports patients who cannot stand completely on their own.
- Pelvic Stabilizer – supports and stabilizes the pelvis.
- Positioning Bar – provides adjustable depth.
- Vertical Handgrips – provides alternate handgrip position, posterior models only.
- Variable Resistance Tabs – controls the speed of the gait trainer and aligns an asymmetrical gait pattern.
- One-Directional Rear Wheels – allows only forward mobility, posterior models only.
- Swivel Wheel Locking Bracket – provides option to lock swivel wheels.
- Limited Lifetime Warranty on the frame.

## ACCESSORIES:

- CE 1035 Set of adjustable forearm platforms, 1pr/bx
- CE 1035 S Set of mounting brackets for forearm platform, 1pr/bx
- CE 1060 Set of ankle prompts, 1pr/bx
- CE 1080 Trunk support for anterior and posterior gait trainers, 1/bx
- CE 1080 S Trunk support for pediatric & tyke anterior & posterior gate trainers, 1/bx
- CE 1080 L Trunk support for adult anterior & posterior gait trainers, 1/bx
- CE 1260 Set of extended uprights for the LT 3100, LT 3200, LT 2100 and LT 2200, 1pr/bx
- CE 1270 Set of extended uprights for the LT 3000 and LT 2000, 1pr/bx

ITEM #	ADJUSTABLE HEIGHT	INSIDE WIDTH BETWEEN HANDLEBAR	WIDTH
LT 3000	29" - 36"	21.5"	26.5"
LT 3100	22" - 28"	16"	23.5"
LT 3200	17" - 20"	14.5"	22"
LT 2000	29" - 36"	21.5"	26.5"
LT 2100	22" - 28"	16"	23.5"
LT 2200	17" - 20"	14.5"	22"